

Omnilux Clear



Clinically proven benefits:

1. Reduces active acne and blemishes, redness & inflammation
2. Helps clear breakouts faster
3. Reduces and helps soften the appearance of acne scarring
4. All natural, no pain, side effects or downtime



Tough on acne. Gentle on skin.

The only flexible LED mask that passes FDA, TGA, and CE regulations for treating acne. Advanced home-use light therapy by Omnilux allows for dermatologist grade treatments right from home, and is the single best step you can take in restoring your acne prone skin and maintaining clear results.



How to use:

Four 10-minute treatments for six weeks. After the initial series is complete, use at the first signs of a breakout.

- Charge your controller & then connect it to the device
- Prepare your skin by cleansing and patting dry
- Secure the device to your head using the velcro straps
- Start the 10-minute treatment
- Remove the mask and apply moisturizer to skin
- Gently wipe down the mask

Average acne clearance after 4 weeks of use:

70%



Real users, real results



^ BEFORE

^ AFTER

^ BEFORE

^ AFTER

^ BEFORE

^ AFTER

Technical specifications:

LED TYPE

Blue and red

WAVELENGTHS

415nm (blue) and 633nm (red)

NUMBER OF LEDS

132 LEDs (2 per LED bulb)

66 LED bulbs

TARGETED USE

Face

CERTIFICATIONS

FDA, CE, TGA

WARRANTY

2 years: device

2 years: controller

INCLUDES

- Flexible silicone Clear device
- Rechargeable battery controller
- Power supply & plug adaptors
- USB C line to USB A connector
- 2x adjustable straps
- Carry bag
- Built-in protective eye shields

Not sure if Omnilux Clear is right for you? Consult with your provider. For more information, please visit

www.omniluxLED.com

+1-707-309-7024

support@omniluxLED.com

